

LYNX HIV p24 Antigen Test

Northwestern Global Health Foundation

Kara Palamountain

LYNX HIV p24 Antigen Test
CE Marked
Lot #
Exp. 2014-12
MICRO
SQUID

LYNX HIV p24 Antigen Test Summary

- Type of assay, principal of test
 - Qualitative p24 antigen based immunochromatographic assay
- Target of test (subtypes)
 - HIV-1
- Specimen types
 - 80µl whole blood
- Reference test
 - Roche Diagnostics COBAS AmpliPrep/COBAS Taqman HIV-1 Test (TNA PCR) alternate approved test (e.g., Abbott)
- Components
 - Instrument, blood collection tube, plasma separator, buffer, test strip, test strip reader

Target setting


- Target patient
 - 4 weeks to 18 months, although studies in infants less than 4 weeks are planned.
- Level of HCC targeted
 - Expected target: DBS collection sites, but studies are planned to assess this
- Type of health care worker suitable for
 - Expected target: All staff levels, but studies are planned to assess this
- Footprint (dimensions and weight)
 - Dimensions: 20.2cm x 15.6cm x 13.4cm / Weight: 1.7kg
- Pluripotency
 - Dedicated instrument used solely for EID; no future tests anticipated
- Cold chain, kit stability and storage
 - No cold chain required; stability studies on locked design still underway
- Training requirements
 - Less than 1 day
- Throughput per day, time to results
 - No batching. The platform can test 11.7 tests per 8 hour day
 - 53 minutes (allowing two minutes for specimen collection)
- Power requirements
 - The platform has a built-in rechargeable battery (up to 8 hours)
- 3rd Party consumables
 - Finger/heel stick consumables (gloves, lancet, alcohol swab, gauze pad)
- Waste
 - Blood collection tube, test strip, plasma separator

Procedure


Step 1: Collect blood




Step 2: Apply blood to LYNX Plasma Separator



Step 3: Plunge Plasma Collection Pad into the Reaction Tube




Step 4: Separate Reaction Tube from LYNX Plasma Separator



10 minutes

Step 5: Add LYNX Buffer



Step 6: Heat



11 minutes

Step 7: Insert LYNX Test Strip



30 minutes

Step 8: Read Test















Product availability

- Cost per test and per equipment, maintenance, ancillary equipment (details of volume tiers)
 - Initial pricing \$2000 and \$500 for reader; lower prices as volumes increase
 - Initial pricing \$15; lower prices as volume increase
- Lock down product
 - As of Nov. 2014, the test format and instrument designs are locked down; a prototype reader is in the final design stage
- Regulatory status and plans
 - CDC and ERPD by 2016; WHO approval by 2017
- Manufacturing capacity (Q1 Q2, Q3, Q4 2015 n= tests and equipment)
 - Will not be available for commercial use by Q1 2015
- Commercial availability?
 - High volume manufacturing in place by Q1 2016 but laboratory and field evaluations still need to be completed
- Global in country support
 - Regional support under development





Technology Performance to Date

- No evaluations have been performed on the latest locked design.
- Prior studies indicate assay performs in laboratory settings
 - Verified on residual specimens being tested by total nucleic acid (TNA) PCR (Roche AmpliPrep/COBAS Taqman HIV-1) at the National Health Laboratory Service Virology Laboratory in Groote Schuur Hospital, Observatory, South Africa.
 - A total of 691 subjects were tested, of which 642 were between the ages of 4 weeks and 18 months. Of these, 80% originated in clinics and 20% in hospitals.
 - TNA PCR detected 32 (5.0%) positives, of which, 30 were positive by the LYNX p24 Test and 2 were equivocal.
 - Classifying equivocals as positive, the sensitivity was 100% (95% confidence interval: 90 - 100%).
 - Classifying equivocals as negative, the sensitivity was 94% (78 - 97%). In both cases, specificity was 99% (99 - 100).
 - Robustness was 100%, since there were no invalids, and turnaround times were less than one hour.
- But user variability in POC settings necessitated product changes

Table 1: Review of Evaluations of LYNX HIV p24 Antigen Test, Northwestern Global Health Foundation

Phase	Phase I	Phase II	Phase IIIa	Phase IIIb
Country, site, n size, # positives	   Evanston, USA n=200, +50 (2015 Q1)	   Cape Town, ZA n=1000, +50 (2015 Q1)	   Cape Town, ZA n=1000, +50 (2015 Q1)	
		 TM Cape Town, ZA + TBD n=1000, +50 (2015 Q4)	 TM Cape Town, ZA + TBD n=3000, +150 (2015 Q4)	 TM TBD n=3000, +150 (2016)

Key

Evaluation with locked down product	TM
Evaluation with Prototype	
Developer led evaluation	
Evaluation In progress	
Evaluation Scheduled (start date)	
No information	NA

LYNX HIV p24 Antigen Test Development Team

Robert Elghanian
Abhishek Agarwal
Zaheer Parpia
Swati Satish
Tanima Jahan Abedin
Lori Johnson-Brand

Kerry Berland
Paul Berland
Mark Kmiecik

Tom Westberg

David Kelso
Mark Fisher

Diana Hardie

